



Resolute

ZOTAROLIMUS-ELUTING CORONARY STENT SYSTEM

RESOLUTE International Trial

12-month results from TCT 2010

Strong New Data from TCT 2010

RESOLUTE International 12-month summary

- **Strength in numbers**
 - RESOLUTE International is a large, high quality all-comer trial
 - 2349 real world patients enrolled in 88 International sites
 - 97.4% patient follow-up at 12 months
 - Excellent clinical outcomes at 12 months
 - Cardiac death/TVMI (Primary Endpoint): 4.1%
 - ARC definite/probable stent thrombosis (Secondary Endpoint): 0.87%
- **Strength in consistency**
 - Total of 3489 real world Resolute DES patients in RESOLUTE All Comers and RESOLUTE International trials
 - Clinical outcomes are consistent between the trials at 12 months
- **Strength in complexity**
 - Resolute DES continues to show strong results in complex patients
 - Almost 70% of patients enrolled were defined as complex
 - Consistent outcomes are seen across patient subgroups

RESOLUTE All Comers and RESOLUTE International were not specifically designed or powered for complex patient subset analysis.

RESOLUTE International Trial Design


Strength in numbers: large, all-comer trial

PI: J. Belardi, F-J. Neumann, P. Widimský

Real World (Open Label)
All patients with symptomatic coronary artery disease eligible for DES implantation (no lesion/vessel limitations)

Resolute Stent
n = 2200

88 international sites (Europe, Asia, Africa and South America)
No angiographic follow-up
100% Independent clinical event adjudication
25% randomly assigned to 100% monitoring



Primary endpoint: Composite of cardiac death and target vessel MI at 12 mo
Key secondary endpoint: ARC definite/probable stent thrombosis at 12 mo
Drug therapy: ASA and clopidogrel/ticlopidine ≥6 mo (per guidelines)



Enrollment Reflects Real-World Practice

Baseline characteristics similar to RESOLUTE All Comers


	<i>RESOLUTE International</i>	<i>RESOLUTE All Comers</i>	
	<i>Resolute DES</i> (n = 2349)	<i>Resolute DES</i> (n = 1140)	<i>Xiience V DES</i> (n = 1152)
Age (yr)	63.5 ± 11.2	64.4 ± 10.9	64.2 ± 10.8
Men (%)	77.8	76.7	77.2
Diabetes mellitus (%)	30.5	23.5	23.4
Insulin Dependent (%)	9	8.4	7.1
Prior MI (%)	27	28.9	30.4
Unstable Angina (%)	26.1	19.4	18.9
AMI (within 12 hr) (%)	9.7	15.4	17.8
AMI (within 72 hr) (%)	20	28.9	28.8
Lesions treated per patient	1.3 ± 0.7	1.5 ± 0.7	1.5 ± 0.8
Multi vessel treated (%)	14.0	25	25
Small vessel (RVD ≤2.75 mm)	45.4	67.8	67.4
Long lesion (length >18 mm)	46.1	18.2	21.2
In-stent restenosis (%)	7.6	8.1	8.0
Bifurcation/trifurcation (%)	18.2	16.9	17.7
Total occlusion (%)	6.3	16.3	17.2
Complex Patients¹ (%)	67.5	67.0	65.6

¹Complex patient definition: bifurcation, bypass graft, ISR, AMI <72 hr, LVEF <30%, unprotected LM, >2 vessels stented, renal insufficiency or failure (creatinine >140 µmol/L), lesion length >27 mm, >1 lesion per vessel, lesion with thrombus or TO (preprocedure TIMI = 0).

With the exception of small vessels and long lesions (treatable with a single 38-mm length stent), Resolute DES currently is not specifically approved for the patient subsets noted above.

Low Cardiac Death/Target Vessel MI Rate

RESOLUTE International primary endpoint




<i>RESOLUTE All Comers</i> 12 months	<i>Resolute DES</i> n = 1119	<i>Xiience V DES</i> n = 1126	<i>p-Value</i>
Cardiac death/TVMI	5.3%	5.4%	0.96

RESOLUTE All Comers was not specifically designed or powered to individually compare cardiac death/TVMI.

Low Stent Thrombosis Rate for an All-Comer Trial


RESOLUTE International secondary endpoint



Stent Thrombosis Rates in All-Comer Trials

Studies not powered for this low frequency ST event


ARC Definite/Probable ST 12-Month Data



For distribution only in markets where Resolute has been approved. Not for distribution in the USA or Japan. © 2010 Medtronic, Inc. All rights reserved. UC2011027478 EE 9/10

Excellent Safety and Efficacy Clinical Outcomes

RESOLUTE International 12-Month Data




Results are Consistent with RESOLUTE All Comers

Strength in consistency: total of 3489 Resolute patients enrolled

Composite Endpoints at 12 Months


- RESOLUTE All Comers (Xience V DES, n = 1126)
- RESOLUTE All Comers (Resolute DES, n = 1119)
- RESOLUTE International (Resolute DES, n = 2287)



TLF: cardiac death, TVMI, clinically driven TLR

RESOLUTE International – A Truly Real-World Trial

Strength in complexity: almost 70% of patients are complex




- Bifurcation
- Bypass grafts
- ISR
- AMI <72 hr
- LVEF <30%
- Unprotected LM
- >2 vessels stented
- Renal insufficiency or failure (creatinine >140 $\mu\text{mol/L}$)
- Lesion length >27 mm
- >1 lesion per vessel
- Lesion with thrombus
- Total occlusion (preprocedure TIMI = 0)

With the exception of long lesions (treatable with a single 38-mm length stent), Resolute DES currently is not specifically approved for the patient subsets noted above.

Resolute DES Shows Strong Results in Complex Patient Subgroups

Complex Patients: Composite Endpoints at 12 Months

- RESOLUTE All Comers (Xience V DES, n = 742)
- RESOLUTE All Comers (Resolute DES, n = 752)
- RESOLUTE International (Resolute DES, n = 1545)




¹Complex patient definition: bifurcation, bypass graft, ISR, AMI <72 hr, LVEF <30%, unprotected LM, >2 vessels stented, renal insufficiency or failure (creatinine >140 µmol/L), lesion length >27 mm, >1 lesion per vessel, lesion with thrombus or TO (preprocedure TIMI = 0). With the exception of long lesions (treatable with a single 38-mm length stent), Resolute DES currently is not specifically approved for the patient subsets noted in this complex patient definition.

RESOLUTE All Comers and RESOLUTE International were not specifically designed or powered for complex patient subset analysis.

Similar Performance Across Complex Patient Subsets

RESOLUTE International 12-Month Data



Resolute DES is not specifically approved for STEMI, multivessel disease and diabetic patient subgroups.

Strong New Data from TCT 2010

RESOLUTE International 12-month summary

- **Strength in numbers**
 - RESOLUTE International is a large, high quality all-comer trial
 - 2349 real world patients enrolled in 88 International sites
 - 97.4% patient follow-up at 12 months
 - Excellent clinical outcomes at 12 months
 - Cardiac death/TVMI (Primary Endpoint): 4.1%
 - ARC definite/probable stent thrombosis (Secondary Endpoint): 0.87%
- **Strength in consistency**
 - Total of 3489 real world Resolute DES patients in RESOLUTE All Comers and RESOLUTE International trials
 - Clinical outcomes are consistent between the trials at 12 months
- **Strength in complexity**
 - Resolute DES continues to show strong results in complex patients
 - Almost 70% of patients enrolled were defined as complex
 - Consistent outcomes are seen across patient subgroups

RESOLUTE All Comers and RESOLUTE International were not specifically designed or powered for complex patient subset analysis.



Resolute

ZOTAROLIMUS-ELUTING CORONARY STENT SYSTEM

A Strong Choice For Your Complex Daily Practice

Medtronic Vascular
3576 Unocal Place
Santa Rosa, CA 95403 USA
Tel: +1.707.525.0111

www.Medtronic.com

Medtronic BV
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel: +31.45.566.8000
Fax: +31.45.566.8668